

Mini-Regulator for natural gas and refinery gas standards

- 0–300 psig inlet pressure range.
- 0–15 psig outlet pressure range.
- Supplied with 0–15 psig outlet pressure gauge, brass CGA 170 nut and nipple, and 1/8-inch compression outlet.

Description	qty.	cat.#
Mini-Regulator	ea.	22032



High-Purity VOC Regulators

- Single-stage, stainless steel.
- Two pressure gauges and CGA-180 fitting.
- Stainless steel diaphragm and Kel-F® seat.
- Accurate pressure control even at low flow rates.
- Individually tested for leaks and impurities.

Spectra Gas 7621 High-Purity VOC Regulators

Specifications: Outlet Connection: 1/8" tube compression
Maximum Inlet Pressure: 3,000 psig Internal Volume: 3.03 cc

Description	qty.	cat.#
0–30 psig outlet pressure gauge	ea.	21572
0–100 psig outlet pressure gauge	ea.	21572-R100



Air Liquide High-Purity VOC Regulators

Specifications: Outlet Connection: 1/4" NPT female
Maximum Inlet Pressure: 3,000 psig Internal Volume: 3.03 cc

Description	qty.	cat.#
CGA-180 (0–30 psig)	ea.	26371
CGA-180 (0–100 psig)	ea.	26372



Regulators for use with 14 L and 48 L Scott/Air Liquide transportable gases

Use the CGA-160 inlet connection with 14 L Scott/Air Liquide transportable gases.
Use the CGA-165 inlet connection with 48 L Scott/Air Liquide transportable gases.

Specifications: Maximum Inlet Pressure: 300 psig
Outlet Pressure Range: 2–10 psig
Maximum Delivery Pressure: 25 psig
Operating Temperature Range: 35 °F to 150 °F (2 °C to 65 °C)
Outlet Connection: 1/4" female NPT

Materials of Construction:
Body: Brass
Diaphragm: Viton®
Seat: Acetal
Seal: Viton®

Description	qty.	cat.#
Regulator, CGA-160 Inlet Connection	ea.	22690
Regulator, CGA-165 Inlet Connection	ea.	22691



also available

Single-Stage and Dual-Stage
Ultra-High Purity Gas Regulators
See **pages 301–303**.



Syringe Adaptor Kit for Single-Stage VOC Regulator

Use to withdraw sample from a high-pressure cylinder after pressure reduction through the high-purity VOC single-stage regulator.

Kit contains one nickel-plated brass 1/4" NPT to female luer fitting, which can be used with an A-2 Luer syringe (cat.# 20162 or 20163), and one stainless steel 1/4" NPT x 1/8" compression fitting with septum (can be used with any syringe needle).

Description	qty.	cat.#
Syringe Adaptor Kit	kit	21118

