



24530

24843

Autosampler Syringes for Shimadzu GCs

- Hamilton and SGE syringes are designed and tested to meet critical autosampler performance needs.
- Needle point styles are designed to withstand multiple, rapid injections through a septum.

Hamilton
for Shimadzu AOC 9 GCs

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	cat.#	qty.	Restek cat.#
5 µL	RN	26s	2"/51 mm	2	75RN	87930	ea.	24617
5 µL	N	26s	2"/51 mm	2	75N	87900	ea.	24938
10 µL	RN	26s	2"/51 mm	2	701RN	80330	ea.	24530

SGE
for Shimadzu AOC 14, 17, 20, and 20i GCs

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	cat.#	qty.	Restek cat.#
0.5 µL	R	23	42 mm	Cone	0.5BR-S-0.63	000445	ea.	24841
0.5 µL	R	26	42 mm	Cone	0.5BR-S-0.47	000440	ea.	24840
5 µL	F	23	42 mm	Cone	5F-S-0.63	001988	ea.	24843
10 µL	R	23	42 mm	Cone	10R-S-0.63	002898	ea.	24845
10 µL*	R	23	42 mm	Cone	10R-S-GT-0.63	002902	ea.	24846
10 µL	R	26	42 mm	Cone	10R-S-0.47	002897	ea.	24844

*Gas-tight syringe.



24749

Replacement Needles for Shimadzu Autosampler Syringes

SGE

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	qty.	Restek cat.#
0.5 µL	23	42 mm	Cone	033745	ea.	24747
10 µL	23	42 mm	Cone	037747	2-pk.	24749
10 µL	26	42 mm	Cone	037745	2-pk.	24748

Guide to Needle Termination Codes

Hamilton:

- (N) Cemented Needle
- (RN) Removable Needle
- (ASN) Autosampler Cemented Needle
- (ASRN) Autosampler Removable Needle
- (TLL) PTFE Luer Lock
- (KH) Knurled Hub
- (LT) Luer Tip
- (LTN) Luer Tip Cemented Needle

SGE:

- (F) Fixed Needle
- (R) Removable Needle
- (LL) Luer Lock
- (LT) Luer Tip