

A-2 Luer Gas-Tight Syringes*

- For use with VOA-CYL 4 standards.
- Push-style on-off button valve.

Use to withdraw a gas sample from a high-pressure cylinder through a high-purity VOC single-stage regulator after attaching the appropriate adaptor (cat.# 21118) to the regulator. Use with luer syringe needles.

Description	qty.	cat.#
100 µL Luer Syringe	ea.	20162
500 µL Luer Syringe	ea.	20163

*Restek's A-2 syringes have a luer tip.



20162

High-purity VOC single-stage regulator is described on page 453 of the Air Sampling Products section.

Replacement Needles for A-2 Luer Gas-Tight Syringes

Description	Needle Gauge	Size	qty.	cat.#
Luer Syringe Needles for 100 µL A2 Gas-Tight Syringes	22	0.028" OD x 0.006" ID x 2"	3-pk.	22783
Luer Syringe Needles for 500 µL A2 Gas-Tight Syringes	22	0.028" OD x 0.012" ID x 2"	3-pk.	22784

Heavy Duty Purge & Trap Syringe Dynatech Precision Sampling

- Heavy-duty glass barrel with metal front and rear flanges.
- PTFE luer lock tip.
- Can fill and empty sparge tubes.
- Sample-Lok valve accepts standard luer lock needles.
- Suitable for purge & trap applications.

Syringe	5 mL cat.#	10 mL cat.#
without Sample-Lok	21206	21209
with Sample-Lok	21208	21207



21208

Luer Lock Valve Adaptors

Use a luer lock valve adaptor to connect a luer lock fitting to a Rheodyne® or Valco® valve.

Description	SGE cat. #	qty.	Restek cat.#
Luer Lock Valve Adaptors	200010	2-pk.	21283



21283

Jumbo Syringe

Clear acrylic syringes, ideal for holding and dispensing large volumes of gas. An adjustable plunger on the O-ring ensures that the syringe is gas-tight over a long period of time. The central port is supplied with a luer lock fitting; the secondary port is supplied with a septum nut. This enables access to the gas sample for adding standards or removing a subsample. The plunger stem is detachable, making sample storage easy.

Volume	Model	SGE cat.#	qty.	Restek cat.#
500 ml	500MAR-LL-GT	009910	ea.	21275
1000 ml	1000MAR-LL-GT	009920	ea.	21276
2000 ml	2000MAR-LL-GT	009930	ea.	21277



21276

Syringe O-Rings

Syringe Volume	SGE cat.#	qty.	Restek cat.#
500 mL	032527	ea.	21278
1,000 mL	032532	ea.	21279



21279

21278

Septa also available (cat.# 27127, page 223).