



24532



24703

Manual Microliter Syringes

- Economical.
- Available with cemented needles (N/F) or removable needles (RN/R).
- Each syringe plunger and barrel assembly is manufactured as one working unit; components are not interchangeable or individually replaceable.

Hamilton 700 Series

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	Hamilton cat.#	Restek qty.	Restek cat.#
5 µL	N	26s	2"/51 mm	2	75N	87900	ea.	24938
5 µL	RN	26s	2"/51 mm	2	75RN	87930	ea.	24617
10 µL	N	26s	2"/51 mm	2	701N	80300	ea.	20174
10 µL	N	26s	2"/51 mm	2	701N	80366	6-pk.	20175
10 µL	RN	26s	2"/51 mm	2	701RN	80330	ea.	24530
25 µL	N	22s	2"/51 mm	2	702N	80400	ea.	24531
10 µL	N	26s	2"/51 mm	5	701NPT5	80339	ea.	24967
25 µL*	RN	22s	2"/51 mm	2	702RN	80430	ea.	24532
50 µL	N	22s	2"/51 mm	2	705N	80500	ea.	24533
50 µL*	RN	22s	2"/51 mm	2	705RN	80530	ea.	24534
100 µL	N	22s	2"/51 mm	2	710N	80600	ea.	24535
100 µL*	RN	22s	2"/51 mm	2	710RN	80630	ea.	24536
250 µL	N	22s	2"/51 mm	2	725N	80700	ea.	24537
250 µL*	RN	22s	2"/51 mm	2	725RN	80730	ea.	24538
500 µL	N	22	2"/51 mm	2	750N	80800	ea.	24539
500 µL*	RN	22	2"/51 mm	2	750RN	80830	ea.	24540

*Replacement RN/R needles available. See page 390.

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	SGE cat.#	Restek qty.	Restek cat.#
5 µL	F	26	50 mm	2	5F**	001000	ea.	24700
5 µL*	R	26	50 mm	2	5R**	001050	ea.	24701
10 µL	F	26	50 mm	2	10F**	002000	ea.	24702
10 µL	F	26	50 mm	2	SK-10F**	002030	6-pk.	24715
10 µL*	R	26	50 mm	2	10R**	002050	ea.	24703
10 µL*	R	26	50 mm	2	SK-10R**	002080	6-pk.	24716
25 µL	F	25	50 mm	2	25F	003000	ea.	24704
25 µL*	R	25	50 mm	2	25R	003050	ea.	24705
50 µL	F	25	50 mm	2	50F	004000	ea.	24706
50 µL*	R	25	50 mm	2	50R	004050	ea.	24707
100 µL	F	25	50 mm	2	100F	005000	ea.	24708
100 µL*	R	25	50 mm	2	100R	005050	ea.	24709
250 µL	F	25	50 mm	2	250F	006000	ea.	24710
250 µL*	R	25	50 mm	2	250R	006050	ea.	24711
500 µL	F	25	50 mm	2	500F	007000	ea.	24712
500 µL*	R	25	50 mm	2	500R	007050	ea.	24713

*Replacement RN/R needles available. See page 390.

**With plunger protection.

Guide to Needle Termination Codes

Hamilton:

- (N) Cemented Needle
- (RN) Removable Needle
- (ASN) Autosampler Cemented Needle
- (ASRN) Autosampler Removable Needle
- (TLL) PTFE Luer Lock
- (KH) Knurled Hub
- (LT) Luer Tip
- (LTN) Luer Tip Cemented Needle

SGE:

- (F) Fixed Needle
- (R) Removable Needle
- (LL) Luer Lock
- (LT) Luer Tip

SuperfleX Flexible Plunger Microliter Syringes

- The SuperfleX syringe features a super-elastic plunger that returns to its original shape if twisted or bent.

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		qty.	Restek cat.#
					Model	cat.#		
5 µL	F	26	50 mm	2	5FX	001100	ea.	24720
10 µL	F	26	50 mm	2	10FX	002100	ea.	24722
10 µL	F	26	50 mm	2	SK-10FX	002130	6-pk.	24718
10 µL	R	26	50 mm	2	SK-10RX	002180	6-pk.	24719

Reinforced Plunger Microliter Syringes

- Reinforcing extends the life of fragile, small diameter plungers.
- Economical because of longer life—even in rugged applications.
- Hand-lapped plunger ensures a perfect seal with the barrel.

Hamilton

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		qty.	Restek cat.#
					Model	cat.#		
5 µL	N	26s	2"/51 mm	2	95	87920	ea.	24541
5 µL*	RN	26s	2"/51 mm	2	95	87925	ea.	24542
10 µL	N	26s	2"/51 mm	2	901	80360	ea.	24543
10 µL*	RN	26s	2"/51 mm	2	901	80370	ea.	24544

*Replacement RN/R needles available. See page 390.

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		qty.	Restek cat.#
					Model	cat.#		
5 µL	F	26	50 mm	2	5F-GP	001400	ea.	24769
5 µL*	R	26	50 mm	2	5R-GP	001450	ea.	24766
10 µL	F	26	50 mm	2	10F-GP	002400	ea.	24767
10 µL*	R	26	50 mm	2	10R-GP	002450	ea.	24768

*Replacement RN/R needles available. See page 390.

