



21339

### Flexible Stainless Steel Hoses

Description	Length	Fittings	qty.	cat. #
Flexible Stainless Steel Hose	36"	1/4" Female NPT	ea.	21339
Flexible Stainless Steel Hose	18"	1/4" Female NPT	ea.	21340



21334

### Hydrogen Flashback Arrestor

- Valve cuts off gas supply in the event of a flashback.
- Can be reset and reused after flashback has occurred.
- Meets safety requirements of BS6158, EN730, ISO 5175, AS 4603 and are UL Listed 23Y5.

#### Specifications:

Body and Internals:	Brass with brass internals
Temperature Sensor:	Polyethylene
Valve spring:	316 stainless steel
Connections:	1/4" female NPT
Maximum pressure:	Hydrogen, 50 psig (3.5 bar)

Description	Fittings	qty.	cat. #
Hydrogen Flashback Arrestor, Brass Body	1/4" Female NPT	ea.	21334



21336

### CGA Fittings

CGA-specified nuts and nipples with internal frit, 1/4-inch NPT nickel-plated brass.

Description	qty.	cat. #
CGA 580 Fitting (N <sub>2</sub> , He, Ar)	ea.	21336
CGA 350 Fitting (Hz, P <sub>2</sub> )	ea.	21337
CGA 590 Fitting (Air)	ea.	21338



23134



23179

### Swagelok® Male Connector, Pipe-to-Tube Fittings

Fitting Type	Size (inches)	Similar to Swagelok	Brass		Stainless Steel	
			qty.	cat. #	qty.	cat. #
Male Connector	1/4" to 1/2" NPT	400-1-4	10-pk.	23134	2-pk.	23184
Male Connector	1/8" to 1/4" NPT	200-1-4	10-pk.	23136	2-pk.	23186
Tube End Reducer	1/4" to 1/8"	200-R-4	5-pk.	23129	2-pk.	23179