



22067

**Viton® O-Rings** for Apex Liners

Description	Max. temp.	qty.	cat.#
Viton O-Rings for APEX liners	250 °C	25-pk.	22067



22354



22355

**Accessories** for ATAS LINEX DMI System

Description	qty.	cat.#
Needle Guide for ATAS LINEX DMI Liner	5-pk.	22354
Microvial/Insert for ATAS LINEX DMI Liner	100-pk.	22355



**Sky**<sup>®</sup>  
Inlet Liners

**True Blue Performance**

Exceptionally inert Sky<sup>®</sup> inlet liners with state-of-the-art deactivation improve trace-level analysis—**and now come with a 100% satisfaction guarantee!\***



See **pages 193-202** or visit [www.restek.com/sky](http://www.restek.com/sky)

\* For details on our 100% satisfaction guarantee, visit [www.restek.com/sky](http://www.restek.com/sky)