



**Injector & Autosampler Supplies** for Thermo HPLC Systems

Description	Model #	Similar to Thermo Part #	qty.	cat.#
Rotor Seal Assembly, Rheodyne 7010	TSP AS100, 300, 1000, 3000, 3500, 8875, and 8880 Autosamplers	1535-4048, 7010-039, 2508-0360	ea.	25481
Syringe Assembly, 250 µL	1000, CM3000, CM3500, AS100	A3588-020	ea.	25482
Syringe, 500 µL	8800 Series	3301-0100	ea.	25483



**Lamps** for Thermo HPLC Systems

Description	Model #	Similar to Thermo Part #	qty.	cat.#
Lamp, UV	Spectrachrom 200, UV100, UV1000, UV150, UV200, UV2000, UV3000	9551-0023	ea.	25484



**Pump Supplies** for Thermo HPLC Systems

Description	Model #	Similar to Thermo Part #	qty.	cat.#
Inlet Check Valve Assembly	8800, 8810, ISOCHROM, P1000, P2000, P4000	A3495-010	ea.	25474
Outlet Check Valve Assembly	8700, 8800, 881, ISOCHROM	A3490-010	ea.	25475
Transducer Check Valve Assembly	P1000, P2000, P4000	A3990-010	ea.	25479
Check Valve Cartridge	3000, 3200, 3500, 4000, 4100, Constametric I, II, and III	900946	ea.	25485
Piston	8800, 8810, ISOCHROM, P1000, P2000, P4000	A3102-010	ea.	25476
Sapphire Plunger	CM3000, CM3200, CM3500, CM4000	801306	ea.	25486
Back-Up Seal	8800, 8810, ISOCHROM, P1000, P2000, P4000	A2963-010	ea.	25477
Plunger Seal, Gold Superseal**	8800, 8810, ISOCHROM, P1000, P2000, P4000	A2962-010	ea.	25478
Plunger Seal Kit, Gold** Includes: SS plunger seal spacer, Rulon plunger seal spacer (2), O-ring, gold plunger seal	LDC Constametric Pumps	801892001	kit	25487
Plunger Seal, Black*	CM3000, CM3200, CM3500, CM4000	206129001	ea.	25488
Plunger Seal, Gold**	CM3000, CM3200, CM3500, CM4000	206156001	ea.	25489
Kel-F Washer	P1000, P2000, P4000	A4848-010	ea.	25480



\*Graphite-filled PTFE; best for organic solvents.

\*\*Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion; best for buffers.